Nordisk Nyhedsbrev

Danmark - Norge - Sverige - Finland

Oktober 2011

Välkommen

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
Artikler i dette nummer:

1st Automatic

- a new model for improved identification of blood gas tests

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Pic. 2. Töölö Hospital

Pic. 3. Töölö Hospital

1. Blood sample is aspirated from arterial catheter

2. safePICO Aspirator's unique barcode is identified with a barcode reader

3. Patient ID is linked to sample ID. Sample is then ready for analysis.

1st Automatic - a new model for improved identification of blood gas tests

By: Elina Ristola, Jaana Kotila and Petra Ylikukkonen
Helsinki University Central Hospital - Helsinki, Finland

The purpose of the study was to improve co-operation between the neurosurgical ICU and the laboratory in Töölö Hospital and to further increase patient safety. Also, the aim was to set up a new working model that would ease the work of nursing staff in the ward and in the lab.

Background

The neurosurgical ICU in Töölö Hospital is a 16-bed ward operating in three separate rooms. Blood gases are measured at least four times a day on each patient with the Radiometer ABL825 analyzer placed in the ward. With the current working model there have been occasional situations where blood gas syringes have been linked to wrong patients. This accidental mix up has led to situations where blood gas results have been sent to wrong patients. In order to avoid such issues, ICU decided, in co-operation with the lab, to pilot as the first ward in Finland Radiom-

eter's 1st Automatic system designed for improving the blood gas analysis processes and for minimizing human errors in blood gas measurements. In big and busy wards the probability of errors naturally increases and the ICU wanted to minimize human errors that may have significant effect on patient care and treatment decisions.

Methods

In order to enable implementation of the new working model, each patient bed was equipped with a barcode reader, and the web-based FlexLink application was

installed on each laptop located next to the patient beds. With the FlexLink software, sample and patient ID's are linked together immediately after sampling. After linking patient and sample ID's, the syringe is brought to the ABL825 analyzer and loaded to one of its three automatic sample loading positions. There is no risk of mixing samples, because each safePICO syringe has a unique barcode that is automatically linked to the correct patient ID by the FlexLink application. The ABL825 analyzer automatically performs the blood gas analysis and transfers results to the patient and the hospital information sys-

tem. This new model saves a lot of time, because nurses do not need go back to the analyzer any more to wait for the tests to be ready and to link results to correct patients.

Results

The neurosurgical ICU arranged a questionnaire concerning the new working model at the ward in February 2011. According to the results, major improvement has been made. Human errors related to patient identification have decreased significantly, or have basically been eliminated completely, which has led to increased patient safety. The new 1st Automatic model has also saved working time, because nurses can now wait for the results at bedside without going to the analyzer. Eliminating false patient identifications has also had a major impact in the routines of the laboratory, because earlier false patient result removals and corrections took quite a lot of working time of the lab staff.

Conclusions

ICU and the laboratory are very satisfied with the new working model based on the use of 1st Automatic, and implementation and use of the system have been very easy. The new model has minimized

patient identification errors and has saved time both at the ward and in the laboratory.

Improvement proposals

A few development proposals were also pointed out in the questionnaire: It is important that the nurses remember to remove patient ID stickers from the bed as soon as the patient has left ICU. Also, attention needs to be paid to positioning the barcode readers at each patient bed in order to allow easy identification process. This could be achieved by installing proper holders for the barcode readers next to each bed.

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Lactate and lactic acidosis

Author: Chris Higgins - Biomedical writer

<http://acutecaretesting.org/lactate-lactic>

Abstract

The integrity and function of all cells depend on an adequate supply of oxygen. Severe acute illness is frequently associated with inadequate tissue perfusion and/or reduced amount of oxygen in blood (hypoxemia) leading to tissue hypoxia. If not reversed, tissue hypoxia can rapidly progress to multiorgan failure and death. For this reason a major imperative of critical care is to monitor tissue oxygenation so that timely intervention directed at restoring an adequate supply of oxygen can be implemented. Measurement of blood lactate concentration has traditionally been used to monitor tissue

oxygenation, a utility based on the wisdom gleaned over 50 years ago that cells deprived of adequate oxygen produce excessive quantities of lactate. The real-time monitoring of blood lactate concentration necessary in a critical care setting was only made possible by the development of electrode-based lactate biosensors around a decade ago. These biosensors are now incorporated into modern blood gas analyzers and other point-of-care analytical instruments, allowing lactate measurement by non-laboratory staff on a drop (100 µL) of blood within a minute or two. Whilst blood lactate concentration

is invariably raised in those with significant tissue hypoxia, it can also be raised in a number of conditions not associated with tissue hypoxia. Very often patients with raised blood lactate concentration (hyperlactatemia) also have a reduced blood pH (acidosis). The combination of hyperlactatemia and acidosis is called lactic acidosis. This is the most common cause of metabolic acidosis. The focus of this article is the causes and clinical significance of hyperlactatemia and lactic acidosis. The article begins with a brief overview of normal lactate metabolism.

Troponin Testing and the Tyranny of Distance

Author: Steven L. Weier - Masters of Science

<http://acutecaretesting.org/Troponin-Testing>

Abstract

Chest pain is a frequent cause of hospitalization and comprises a significant proportion of the acute medical workload in the western world. This issue is compounded for those living in regional or remote communities due to limited access

to pathology services and tertiary-referral institutions.

Evidence suggests that the replacement of traditional laboratory-based testing for cardiac troponin with a point-of-care test is associated with improved clinical man

agement of patients presenting with chest pain and leads to a significant and sizeable reduction in duration of stay for these patients. Minimizing the in-hospital stay of patients with chest pain, within safe limits, is crucial in reducing the cost of health care.

Automating blood gas improves preanalytical quality

Authors: Agnes Ivanov - quality manager, United Laboratories of Tartu University Hospital
Anu Tamm - director, United Laboratories of Tartu University Hospital

<http://acutecaretesting.org/preanalytical-quality>

Abstract

Tartu University Hospital has improved the quality of blood gas testing in its clinical laboratory and critical care units by introducing an automatic system for blood gas analysis [1].

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